DE - Pedaltrainer faltbar:

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EN - Pedal Exerciser foldable:

Effective home training for comfortable use while seated. Easy and joint-friendly training for arms and legs.

Safety instructions:

This article is not a toy. Keep it away from children and pets. Place this product on a flat, level and firm floor. Make sure the area is clear to avoid injury. It is recommended to place a carpet or mat under the device to protect the floor. Do not wear loose clothing or jewelry when using the product since it might get caught into the moving parts. If you experience pain, dizziness, or discomfort, stop using the device immediately. Stop using the product if there are any signs of damage or defects. Do not stand on the mini bike or apply full body weight to the exerciser.

- 1 Assembly instructions:** Screw the rotary knob onto the training device as shown in Figure 1. To adjust the resistance individually, loosen or tighten the rotary knob.
- 2 Changing the battery:** To access the batteries, you will need a flat tool to open the back of the digital display panel. Slide out the battery and replace it with an LR44 button cell battery. See figure 2. (Note: Do not pull the control panel too far away from the pedal exerciser or the cable may break.)
- 3 Digital display:** When using the device, press the function button in sequence to scroll through the display options. See figure 3.

To fold up the pedal trainer, pull on the inner rings and then fold the support legs inwards.

Note: It is recommended that the angle between the thigh and calf is less than 120°. Please sit on a chair of the correct height when exercising.

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